


User Instructions for Hospital Bed Side Rails

Side rails DO NOT fall within any weight limitations. Side rails can be deformed or broken if excessive side pressure is exerted on the side rails by pulling or leaning on them. Side rails are for the purpose of preventing an individual from accidentally rolling out of the hospital bed. Side rails are not intended to be used for restraint purposes or pulling oneself out of bed. Misuse can lead to damage to the rail and harm to the client.



1. Loosen "A" almost fully.
2. Insert "B" over top of bed frame.
3. Insert "C" on underside of the top portion of the bed frame.
4. Tighten "A" fully. Bed rail is now ready for use.
5. To lower bed rail, raise the rail slightly with one hand and at the same time pull outward on section "D" and release the rail.
6. To raise rail, lift the rail until section "D" drops under bottom rail. You may have to move "D" out of the way of the bottom of the rail first.