




User Instructions for a Flimm-Fighter Percussor




The applicator handle of your Flimm-Fighter is designed for the shape and size of your hand. To self-apply, hold the applicator handle in your hand and turn the machine on. (Do not turn the machine on unless you are holding the applicator handle. Otherwise, the applicator handle will vibrate without restraint.)

With the applicator handle in your hand; turn the machine on, then place the applicator handle head on your chest or side. Be sure to be in proper postural drainage position prior to applying the self-applicator. For proper postural drainage positions, Flimm-Fighter duration and self-care please consult your health care professional.


Care and Maintenance

Your Flimm-Fighter required no user mechanical maintenance or lubrication ever. If your Flimm-Fighter needs servicing, please call MDA for assistance.


The plastic motor housing and the plastic applicator handle can be cleaned with any normal household spray cleanser. However, never immerse the motor housing or the applicator head or handle in any liquid.




Please do **not** tie the flexible drive cable in a knot.




Please do **not** crimp or sharply bend the flexible drive cable.




Please do **not** forcibly pull on the flexible drive cable.



Please do **not** forcibly pull on the electric power cord.



The motor housing should always face the direction of the patient so as to avoid extreme, sharp bends in the flexible drive cable. Please do **not** operate your Flimm-Fighter with any sharp bend in the flexible drive cord.



Never immerse the motor housing or the applicator handle in any liquid!